

---

## New Product Lines

---

### Catheter Extension Sets with InVision-Plus® Neutral® I.V. Connector System

---

NOTE: All catheter extension sets are made of soft, kink-resistant, NON-DEHP PVC Tubing

#### 7" Macrobore Catheter Extension Set with pre-attached InVision-Plus Neutral I.V. Connector

Product Profile

Reorder Number


**InVision-Plus Neutral I.V. Connector**

**RYM-5307**

Slide clamp, rotating collar male-luer lock.

Note: I.V. connector not bonded to set

Tubing inside diameter equivalent to a 16g catheter flow rate



Packaging: 100 Units per Box – 4 Boxes per Case – TOTAL: 400 UNITS PER CASE

---