




Intraluminal Protection System InVision-Plus[®] RED[™] with Neutral Advantage[™] Technology



For Visual Identification of Arterial Catheters or Non-Venous Catheters

- Use with arterial catheters and stopcocks
- Addresses “tubing misconnections - a persistent and potentially deadly occurrence”

The Joint Commission Sentinel Event, Issue 36, April 3, 2006



Fail-safe engineering design supports successful outcomes

	10-Point Standard	InVision-Plus [®]
Infection Barriers	1 Smooth, swabbable septum surface	✓
	2 Septum seal integrity (with no gaps or openings)	✓
	3 Patented double microbial barrier	✓
Outcomes Enhancements	4 Straight-through fluid pathway	✓
	5 Zero dead space	✓
	6 Zero fluid displacement	✓
	7 Low priming volume at 0.027mL	✓
Clinical Benefits	8 100% effective blood clearing at 4.0mL	✓
	9 Saline-only flushing option (Heparin-free)	✓
	10 No clamping sequence or positive pressure syringe technique required	✓



Catalog Number RYM-2001: Please see the RyMed Product Catalog at www.rymedtech.com for the entire family of RyMed's products with Neutral Advantage[™] Technology